

Haier Biomedical

Intelligent Protection of Life Science

-40°C Biomedical Freezer

DW-40W100



Ergonomic Design



Lockable design



Intelligent alarm



Test hole design



Three Alarm Functions

High/low temperature, power failure, and sensor error



Adjustable Temperature

Adjustable freezer compartment range: -20 °C ~ -40 °C

Specifications

DW-40W100

Technical Data

Type	-40°C Deep Freezer
Cabinet Type	Chest
Climate Class	4
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB(A))	≦ 45

Performance

Cooling performance(°C)	-40
Temperature Range(°C)	-20~-40

Control

Controller	Microprocessor
Display	LCD

Electrical Data

Power Supply(V/Hz)	220-240/50	115/60
Power(W)	260	255
Electrical Current(A)	2.5	5.0

Dimensions

Capacity(L/Cu.Ft)	100/3.5
Net/Gross Weight(approx)	43/46(kg) 94.8/101.4(lbs)
Interior Dimensions(W*D*H)	500*385*630(mm) 19.7*15.2*24.8(in)
Exterior Dimensions(W*D*H)	685*650*810(mm) 27.0*25.6*31.9(in)
Packing Dimensions(W*D*H)	740*650*880(mm) 29.1*25.6*34.6(in)
Container load (20'/40'/40'H)	54/108/162



Functions

High/Low Temperature	Y
Power Failure	Y
Sensor Error	Y

Accessories

Test Hole	Y
-----------	---

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R.China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.

