

# Haier Biomedical

Intelligent Protection of Life Science

## -40°C Biomedical Freezer

DW-40W380



### Ergonomic Design



Lockable design



Intelligent alarm



### Three Alarm Functions

High/low temperature, power failure, and sensor error



### Adjustable Temperature

Adjustable freezer compartment range: -20 °C ~ -40 °C

# Specifications

## DW-40W380

### Technical Data

Type	-40°C Deep Freezer
Cabinet Type	Chest
Climate Class	4
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB(A))	≦ 45

### Performance

Cooling performance(°C)	-40
Temperature Range(°C)	-20~-40

### Control

Controller	Microprocessor
Display	LCD

### Electrical Data

Power Supply(V/Hz)	220-240/50	208-230/60
Power(W)	550	550
Electrical Current(A)	5.1	5.1

### Dimensions

Capacity(L/Cu.Ft)	380/13.4
Net/Gross Weight(approx)	82/89kg 180.8/196.2(lbs)
Interior Dimensions(W*D*H)	1375*453*640(mm) 54.1*17.8*25.2(in)
Exterior Dimensions(W*D*H)	1554*724*838(mm) 61.2*28.5*33.0(in)
Packing Dimensions(W*D*H)	1635*760*900(mm) 64.4*29.9*35.4(in)
Container load (20'/40'/40'H)	18/42/63



### Functions

High/Low Temperature	Y
Power Failure	Y
Sensor Error	Y

### Others

Certificate	CE
-------------	----

## Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1  
Haier Road, Qingdao, 266101, P.R.China  
Tel: +86-0532-88936011 / 5955  
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.

