

Haier



MEDICAL INFRARED THERMOMETER

RAPID INDUCTION ACCURATE MEASUREMENT

AET-R1B1

- Highly sensitive infrared probe
- Noncontact measurement
- 1 second can be measured

LED Digital Display

The reading in the black display box is clear and easy to obtain the temperature value




AET-R1B1

One machine multi-purpose

The human body/environment mode is automatically switched which can measure the forehead and the ambient temperature giving more comprehensive care



measure forehead temperature



measure environmental temperature

1 second measurement

High sensitive infrared temperature measurement the temperature value can be measured in one second. It is suitable for multi person temperature measurement or multiple temperature measurement comparison



Ergonomic Design Comfortable to Hold

Following ergonomic design rules the body is solid and stable, handiness comfortable grip, firm grasp



Authoritative Certification of CE/FDA/CFDA

CFDA CE FDA



Features

- One button measurement
- memory storage
- replaceable battery
- forehead measurement
- hot tip
- LED display
- measure in 1s
- °C/°F switch
- intelligent accurate
- automatic shutdown

Technical parameters

Measuring range: 32.0°C~42.2°C
Accuracy: ±0.2°C (35.0°C~42.0°C)
±0.4°F (95.0°F~107.6°F)
±0.3°C (±0.5°F) other range
Resolution: 0.1°C

Memory: 32 groups
Power supply: 2*AAA(DC 3V)
Dimension: 144mm*101mm*26mm
Weight: 88g(exclude battery)
Auto power-off: 60s