

Haier Biomedical Equipment Specifications

Product	ULT freezer model	Rack model	Rack and box discription	Freezer Rack Quantity	Freezer Box Quantity	Cryo Tube Quantity (2ml)
Racks for ULT Freezers	DW-86L338J/DW-86L338JA	DCJ-54-A (SUS304)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	6	120	12000
		DCJ-54-B (SUS201/304) (Middle, 2 layers)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	6	120	12000
		DCJ-44-A (SUS304) (Top and bottom layers)	2 inch stainless steel side-access freezer rack 4*4 (2 inch freezer box*16 units)	6	96	9600
		DCJ-44-B (SUS201/304) (Top and bottom layers)	2 inch stainless steel sliding-drawer freezer rack 4*4 (2 inch freezer box*16 units)	6	96	9600
	DW-86L388J/418S/418ST/419	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	12	300	30000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	12	300	30000
	DW-86L486E	DCJ-54-A (SUS304)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	16	320	32000
		DCJ-54-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	16	320	32000
	DW-86L490J/490JA	DCJ-54-A (SUS304) (Top 2 layers)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	8	160	16000
		DCJ-54-B (SUS201/304) (Top 2 layers)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	8	160	16000
		DCJ-44-A (SUS304) (Bottom 2 layers)	2 inch stainless steel side-access freezer rack 4*4 (2 inch freezer box*16 units)	8	128	12800
		DCJ-44-B (SUS201/304) (Bottom 2 layers)	2 inch stainless steel sliding-drawer freezer rack 4*4 (2 inch freezer box*16 units)	8	128	12800
	DW-86L628E	DCJ-54-A (SUS304)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	20	400	40000
		DCJ-54-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	20	400	40000
	DW-86L578J/578JA/578S/578ST/579/579BP/579BPT	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	16	400	40000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	16	400	40000
	DW-86L728J/728JA/728S/728ST/728SAT/729/729BP/729BPT	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	20	500	50000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	20	500	50000
	DW-86L828J/828JA/828S/828ST/829/829BP/829BPT/828W	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	24	600	60000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	24	600	60000
	DW-86L959BP/959BPT/959W	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	28	700	70000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	28	700	70000
	DW-86L100J	DCJ-53-A (SUS304)	2 inch stainless steel side-access freezer rack 5*3 (2 inch freezer box*15 units)	4	60	6000
		DCJ-53-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*3 (2 inch freezer box*15 units)	4	60	6000
DW-86W420J/420JA	DCJ-10-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	27	270	27000	
DW-86W100J	DCJ-08-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*8 units)	9	72	7200	
DW-150W200	DCJ-10-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	12	120	12000	
DW-150W209	DCJ-10-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	12	120	12000	
	DCJ-12-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*12 units)	12	144	14400	
DW-86L51J	DCJ-53-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*3 (2 inch freezer box*15 units)	2	30	3000	

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
E-mail: inquiry@haierbiomedical.com
Website: www.haiermedical.com



*Haier Biomedical reserves the right to change products and specifications without prior notice.

Haier Biomedical
Providing Superior Global Solutions

